

## Use Cases

UC-01	<i>AED Deployment in Cardiac Emergency Use Case</i>
Description	The base use case for administering an AED.
Primary Actor	AED Operator
Pre-condition	<ol style="list-style-type: none"> <li>1. An AED should only be utilized under the condition that the patient exhibits signs of a suspected cardiac arrest, which include unconsciousness, absence of normal breathing, and lack of a pulse or other signs of circulation. Furthermore, the device must be operated by individuals who have received training in basic life support/AED, advanced life support, or a similar physician-authorized emergency medical response program.</li> <li>2. The device is off</li> </ol>
Post-condition	The AED performs all intended actions properly.
Main Sequence	<ol style="list-style-type: none"> <li>1. Turn the AED unit on by pressing the power button. The AED unit performs an automatic self-test to ensure that the battery usage indicator shows adequate battery capacity for usage and that the defibrillation electrodes are properly connected.</li> <li>2. Check to see if the self-test unit goes from an "X" to a checkmark and the AED will say "Automatic Defibrillator Unit OK".</li> <li>3. AED will output a message saying, "Stay calm".</li> <li>4. The indicator light on the "check responsiveness" picture flashes.</li> <li>5. AED outputs a message saying, "Check responsiveness".</li> <li>6. AED operator Shouts, "Are you OK". They should be shaking the person while doing this.</li> <li>7. The indicator light on the "call for emergency" starts flashing.</li> <li>8. AED outputs a message saying, "Call for help".</li> <li>9. AED operator calls or sends someone else to call for help.</li> <li>10. The indicator light on the "attach electrode pads" picture starts flashing.</li> <li>11. AED outputs a message saying, "Attach defib pads to patient's bare chest".</li> <li>12. AED operator attaches adult or infant/child electrode pads to the patient's bare chest. Follow electrode pad package instructions.</li> <li>13. The indicator light on the "analyzing patient" picture starts flashing.</li> <li>14. AED outputs a message saying, "Don't touch the patient. Analyzing".</li> <li>15. The indicator light on the "perform CPR" picture starts flashing.</li> <li>16. AED Administrator administers CPR. Give two breaths for every thirty compressions.</li> <li>17. The AED will output feedback messages based off your compressions.</li> <li>18. Keep repeating steps 14-16.</li> </ol>

Extensions	<p>1a. If at any point emergency services arrive, the use case ends.</p> <p>1b. If at any point the electrodes are not attached the AED will stop where its at and output a message saying, "Attach electrode pads".</p> <p>1c. if at any point the user presses the power off button go to UC-04.</p> <p>2a. The AED self-test unit stays as an "X" because there is not enough battery life, go to UC-02.</p> <p>2b. When the victim is less than 8 years of age or weighs less than 55 lbs (25 kg), the Fully Automatic AED Plus should be used with ZOLL AED Plus Pediatric Electrodes.</p> <p>2c. The AED self-test unit stays as an "X" because the defibrillation electrodes are not properly attached. The AED will output a message saying, "Fully Attach Defibrillation electrodes to the device".</p> <p>14a. AED detects a shockable rhythm. Go to UC-03. This shockable rhythm could either be the patient has ventricular fibrillation or pulseless ventricular tachycardia.</p> <p>14b. AED outputs a message saying, "No shock advised.".</p> <p>14c. If the shock revives the patient and the patient has a normal heartbeat the use case ends.</p> <p>16a. If the shock revives the patient and the patient has a normal heartbeat use case ends.</p>
Variation	<p>16. When giving compressions the AED will output either one of these 3 messages:</p> <ul style="list-style-type: none"> <li>a) The compression is too deep. Reduce force.</li> <li>b) Compression is too shallow. Press harder.</li> <li>c) Good compression depth. Keep going.</li> </ul>

UC-02	<i>AED Battery Change Use Case</i>
Description	This use case describes a situation where the AED does not have enough battery power to turn on, so a battery change is performed.
Primary Actor	AED operator, Battery
Pre-condition	AED does not have enough battery power to turn on.
Post-condition	AED has a full battery.
Main Sequence	<ol style="list-style-type: none"> <li>1. Ensure that the Fully Automatic AED Plus is turned off. Open the battery compartment by removing the battery cover from the back of the unit.</li> <li>2. Insert 10 new batteries into the Fully Automatic AED Plus unit.</li> <li>3. Connect electrode cable to Fully Automatic AED Plus unit and pack sealed. <ol style="list-style-type: none"> <li>a) electrodes inside unit cover. Close cover.</li> </ol> </li> <li>4. Turn the unit on and wait for the "Unit OK" audio message. Verify that the unit issues appropriate "Adult Pads" or "Pediatric Pads" audio messages.</li> <li>5. Turn the unit off.</li> </ol>

	<ol style="list-style-type: none"> <li>6. Wait 2 minutes. Verify that the green check symbol appears in the status indicator window (located on the left side of the handle) and that the unit does not emit a beeping tone.</li> <li>7. Place a Fully Automatic AED Plus unit in service.</li> <li>8. Check the Fully Automatic AED Plus unit periodically to ensure that the green check symbol appears in the status indicator window</li> </ol>
--	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

UC-03	<i>AED Advises Shock Use Case</i>
Description	This use case describes the situation where the AED determines that a shock is advisable.
Primary Actor	AED Operator
Pre-condition	AED operator has completed use case 1 steps 1-11 properly.
Post-condition	AED administers a shock to the patient.
Main Sequence	<ol style="list-style-type: none"> <li>1. The button in the middle starts flashing.</li> <li>2. The AED operator presses the button in the middle which contains a heart.</li> <li>3. AED outputs a message saying, "Don't touch the patient. Analyzing".</li> <li>4. AED Operator gives a warning to other bystanders "Stand clear".</li> <li>5. AED outputs a message saying, "Shock will be delivered".</li> <li>6. AED shock tone beeps and a shock to the patient is delivered.</li> </ol>
Extensions	1a. The AED does not have enough power to administer a shock. Go to UC-02

UC-04	<i>Power Off Use Case</i>
Description	The concrete use case for power off of the device.
Primary Actor	AED Operator, Battery
Pre-condition	AED has been powered on.
Post-condition	AED has been powered off.
Main Sequence	<ol style="list-style-type: none"> <li>1. The user presses the power button.</li> <li>2. The AED machine powers off.</li> </ol>

